

Activation-induced cytidine deaminase

Catalog No: #AP80386



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Activation-induced cytidine deaminase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Cytidine aminohydrolase,AID,Cytidine aminohydrolase
Accession No.	Uniprot ID:Q75W64
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.