

Adiponectin

Catalog No: #AP80398



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adiponectin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Adipocyte-specific protein AdipoQ, Adipocyte, C1q and collagen domain-containing protein, 30 kDa adipocyte complement-related protein, Adipocyte complement-related 30 kDa protein (ACRP30)
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.