

Apolipoprotein L1

Catalog No: #AP80400



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein L1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apolipoprotein L,Apo-L,Apolipoprotein L-I,ApoL-I,APOL
Accession No.	Uniprot ID:O14791
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.