

Aquaporin-1

Catalog No: #AP80402



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aquaporin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AQP-1,Aquaporin-CHIP,Delayed early response protein 2,DER2,Water channel protein for red blood cells and kidney proximal tubule
Accession No.	Uniprot ID:Q02013
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.