

Aquaporin-1

Catalog No: #AP80417



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Aquaporin-1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Bovine |
| Other Names | AQP-1,Aquaporin-CHIP,Water channel protein CHIP29,Water channel protein for red blood cells and kidney proximal tubule |
| Accession No. | Uniprot ID:P47865 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.