

# Ankyrin repeat domain-containing protein 30A

Catalog No: #AP80425



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Ankyrin repeat domain-containing protein 30A           |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Hu   |
| Other Names        | Serologically defined breast cancer antigen NY-BR-1    |
| Accession No.      | Uniprot ID:Q9BXX3                                      |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.