

Adenylate cyclase type 10

Catalog No: #AP80438



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Adenylate cyclase type 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AH-related protein, Adenylate cyclase homolog, Germ cell soluble adenylyl cyclase, hsAC, Testicular soluble adenylyl cyclase, SAC
Accession No.	Uniprot ID:Q96PN6
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.