

Beta-2 adrenergic receptor

Catalog No: #AP80439



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Beta-2 adrenergic receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta-2 adrenoreceptor,Beta-2 adrenoceptor,Adrb2r
Accession No.	Uniprot ID:P10608
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.