

# AP-1 complex subunit mu-1

Catalog No: #AP80454



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AP-1 complex subunit mu-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AP-mu chain family member mu1A,Adaptor protein complex AP-1 subunit mu-1,Adaptor-related protein complex 1 subunit mu-1,Clathrin assembly protein complex 1 mu-1 medium chain 1,Golgi adaptor HA1/AP1 adaptin mu-1 subunit,Mu-adaptin 1,Mu1A-adaptin
Accession No.	Uniprot ID:Q2KJ81
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.