

# AP-2 complex subunit beta

Catalog No: #AP80455



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AP-2 complex subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AP105B,Adaptor protein complex AP-2 subunit beta,Adaptor-related protein complex 2 subunit beta,Beta-2-adaptin,Beta-adaptin,Clathrin assembly protein complex 2 beta large chain,Plasma membrane adaptor HA2/AP2 adaptin beta subunit,Clapb1
Accession No.	Uniprot ID:P62944
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.