

AP-2 complex subunit beta

Catalog No: #AP80466



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	AP-2 complex subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	AP105B,Adaptor protein complex AP-2 subunit beta,Adaptor-related protein complex 2 subunit beta,Beta-2-adaptin,Beta-adaptin,Clathrin assembly protein complex 2 beta large chain,Plasma membrane adaptor HA2/AP2 adaptin beta subunit,CLAPB1
Accession No.	Uniprot ID:P63009
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.