

Allograft inflammatory factor 1

Catalog No: #AP80469



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Allograft inflammatory factor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AIF-1,Ionized calcium-binding adapter molecule 1,MRF-1,Microglia response factor,Iba1,Mrf1
Accession No.	Uniprot ID:P55009
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.