

# Aldehyde dehydrogenase, dimeric NADP-preferring

Catalog No: #AP80483



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aldehyde dehydrogenase, dimeric NADP-preferring
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ALDHIII,Aldehyde dehydrogenase 3,Aldehyde dehydrogenase family 3 member A1,ALDH3
Accession No.	Uniprot ID:P30838
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.