

Methionine aminopeptidase 2

Catalog No: #AP80488



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methionine aminopeptidase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MAP 2, MetAP 2, Initiation factor 2-associated 67 kDa glycoprotein, Peptidase M 2, p67eIF2, p67, MNPEP, P67EIF2
Accession No.	Uniprot ID:P50579
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.