

Protein argonaute-2

Catalog No: #AP80502



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein argonaute-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Argonaute2,Argonaute RISC catalytic component 2,Eukaryotic translation initiation factor 2C 2,eIF-2C 2,Golgi ER protein 95 kDa,GERp95,Protein slicer,Eif2c2
Accession No.	Uniprot ID:Q9QZ81
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.