

Atrial natriuretic peptide receptor 1

Catalog No: #AP80503



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Atrial natriuretic peptide receptor type A,ANP-A,Guanylate cyclase A,GC-A
Accession No.	Uniprot ID:P18910
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.