

DNA-(apurinic or apyrimidinic site) lyase

Catalog No: #AP80505



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DNA-(apurinic or apyrimidinic site) lyase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	APEX nuclease, APEN, Apurinic-apyrimidinic endonuclease 1, AP endonuclease 1, REF-1, Redox factor-1, APE, APE1, APEX, APX, HAP1, REF1
Accession No.	Uniprot ID:P27695
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.