

Angiopoietin-related protein 4

Catalog No: #AP80538



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Angiopoietin-related protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	425O18-1,Angiopoietin-like protein 4,Fasting-induced adipose factor,Hepatic fibrinogen/angiopoietin-related protein,HFARP,Secreted protein Bk89,Farp,Fiaf,Ng27
Accession No.	Uniprot ID:Q9Z1P8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.