

Acetylserotonin O-methyltransferase

Catalog No: #AP80566



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acetylserotonin O-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Hydroxyindole O-methyltransferase,Hiomt
Accession No.	Uniprot ID:B3GSH5
Tag Info	His
Formulation	Storage Buffer:20mM Tris,0.5mM EDTA,0.5M Arginine,20%glycerol (PH8.0)
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.