

E3 ubiquitin-protein ligase AMFR

Catalog No: #AP80573



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	E3 ubiquitin-protein ligase AMFR
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Autocrine motility factor receptor, isoform 2,AMF receptor, isoform 2,RING finger protein 45,gp78,RNF45
Accession No.	Uniprot ID:Q9UKV5
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.