

# Anti-Muellerian hormone type-2 receptor

Catalog No: #AP80577



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Anti-Muellerian hormone type-2 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Anti-Muellerian hormone type II receptor,AMH type II receptor,MIS type II receptor,MISRII
Accession No.	Uniprot ID:Q8K592
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.