

Annexin A3

Catalog No: #AP80584



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Annexin A3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	35-alpha calcimedin,Annexin III,Annexin-3,Inositol 1,2-cyclic phosphate 2-phosphohydrolase,Lipocortin III,Placental anticoagulant protein III,PAP-III,ANX3
Accession No.	Uniprot ID:P12429
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.