

# Type-1 angiotensin II receptor

Catalog No: #AP80604



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Type-1 angiotensin II receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AT1AR,AT1BR,Angiotensin II type-1 receptor,AT1,AGTR1A,AGTR1B,AT2R1,AT2R1B
Accession No.	Uniprot ID:P30556
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.