

Apolipoprotein C-I

Catalog No: #AP80614



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Apolipoprotein C-I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Apo-CI, Apolipoprotein C1, Liver regeneration-related protein LRRG04
Accession No.	Uniprot ID:P19939
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.