

RAC-alpha serine/threonine-protein kinase

Catalog No: #AP80636



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RAC-alpha serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	AKT1 kinase,Protein kinase B,PKB,Protein kinase B alpha,PKB alpha,Proto-oncogene c-Akt,RAC-PK-alpha,Thymoma viral proto-oncogene,Akt,Rac
Accession No.	Uniprot ID:P31750
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.