

RAC-alpha serine/threonine-protein kinase

Catalog No: #AP80637



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RAC-alpha serine/threonine-protein kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Protein kinase B,PKB,Protein kinase B alpha,PKB alpha,RAC-PK-alpha
Accession No.	Uniprot ID:P47196
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.