

Angiotensin-converting enzyme 3

Catalog No: #AP80645



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Angiotensin-converting enzyme 5, ANG-5, Angiotensin-converting enzyme-like protein 3, ANGPT5, UNQ153/PRO179
Accession No.	Uniprot ID:Q9Y5C1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.