

Angiotensin-converting enzyme 2

Catalog No: #AP80649



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Angiotensin-converting enzyme 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Angiotensin-converting enzyme 2, Heparin-binding epidermal growth factor-like protein 2, HEGF-2, HGF-2, HGF-2L, HGF-2L1, HGF-2L2, HGF-2L3, HGF-2L4, HGF-2L5, HGF-2L6, HGF-2L7, HGF-2L8, HGF-2L9, HGF-2L10, HGF-2L11, HGF-2L12, HGF-2L13, HGF-2L14, HGF-2L15, HGF-2L16, HGF-2L17, HGF-2L18, HGF-2L19, HGF-2L20, HGF-2L21, HGF-2L22, HGF-2L23, HGF-2L24, HGF-2L25, HGF-2L26, HGF-2L27, HGF-2L28, HGF-2L29, HGF-2L30, HGF-2L31, HGF-2L32, HGF-2L33, HGF-2L34, HGF-2L35, HGF-2L36, HGF-2L37, HGF-2L38, HGF-2L39, HGF-2L40, HGF-2L41, HGF-2L42, HGF-2L43, HGF-2L44, HGF-2L45, HGF-2L46, HGF-2L47, HGF-2L48, HGF-2L49, HGF-2L50, HGF-2L51, HGF-2L52, HGF-2L53, HGF-2L54, HGF-2L55, HGF-2L56, HGF-2L57, HGF-2L58, HGF-2L59, HGF-2L60, HGF-2L61, HGF-2L62, HGF-2L63, HGF-2L64, HGF-2L65, HGF-2L66, HGF-2L67, HGF-2L68, HGF-2L69, HGF-2L70, HGF-2L71, HGF-2L72, HGF-2L73, HGF-2L74, HGF-2L75, HGF-2L76, HGF-2L77, HGF-2L78, HGF-2L79, HGF-2L80, HGF-2L81, HGF-2L82, HGF-2L83, HGF-2L84, HGF-2L85, HGF-2L86, HGF-2L87, HGF-2L88, HGF-2L89, HGF-2L90, HGF-2L91, HGF-2L92, HGF-2L93, HGF-2L94, HGF-2L95, HGF-2L96, HGF-2L97, HGF-2L98, HGF-2L99, HGF-2L100
Accession No.	Uniprot ID:Q6TMA8
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.