

Disintegrin and metalloproteinase domain-containing protein 8

Catalog No: #AP80650



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Disintegrin and metalloproteinase domain-containing protein 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ADAM 8,Cell surface antigen MS2,Macrophage cysteine-rich glycoprotein,CD156a,Ms2
Accession No.	Uniprot ID:Q05910
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.