

Serine-protein kinase ATM

Catalog No: #AP80674



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serine-protein kinase ATM
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Ataxia telangiectasia mutated homolog,A-T mutated homolog
Accession No.	Uniprot ID:Q6PQD5
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.