

Atrial natriuretic peptide receptor 2

Catalog No: #AP80700



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Atrial natriuretic peptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Atrial natriuretic peptide receptor type B, ANP-B, Guanylate cyclase B, GC-B
Accession No.	Uniprot ID:P46197
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.