

Aprataxin and PNK-like factor

Catalog No: #AP80701



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Aprataxin and PNK-like factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apurinic-apyrimidinic endonuclease APLF,PNK and APTX-like FHA domain-containing protein,XRCC1-interacting protein 1,C2orf13,PALF,XIP1
Accession No.	Uniprot ID:Q8IW19
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.