

Interferon-inducible protein AIM2

Catalog No: #AP80752



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interferon-inducible protein AIM2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Interferon-inducible protein 210,Ifi-210,Interferon-inducible protein p210,Gm1313,Ifi210
Accession No.	Uniprot ID:Q91VJ1
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.