

# A-kinase anchor protein 17A

Catalog No: #AP80753



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	A-kinase anchor protein 17A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	AKAP-17A,721P,B-lymphocyte antigen,Protein XE7,Protein kinase A-anchoring protein 17A,PRKA17A,Splicing factor, arginine/serine-rich 17A,CXYorf3,DXYS155E,SFRS17A,XE7
Accession No.	Uniprot ID:Q02040
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.