

Activity-regulated cytoskeleton-associated protein

Catalog No: #AP80772



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Activity-regulated cytoskeleton-associated protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	ARC/ARG3.1,Activity-regulated gene 3.1 protein homolog,Arg3.1
Accession No.	Uniprot ID:Q8AWC3
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.