

# Plasma membrane calcium-transporting ATPase 1

Catalog No: #AP80780



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Plasma membrane calcium-transporting ATPase 1          |
| Brief Description  | Recombinant Protein                                    |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE                           |
| Species Reactivity | Pig  |
| Other Names        | PMCA1, Plasma membrane calcium pump isoform 1          |
| Accession No.      | Uniprot ID:P23220                                      |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.