

Plasma membrane calcium-transporting ATPase 2

Catalog No: #AP80789



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Plasma membrane calcium-transporting ATPase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2, PMCA2
Accession No.	Uniprot ID:Q01814
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.