

# Sarcoplasmic/endoplasmic reticulum calcium ATPase

1

Catalog No: #AP80825



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sarcoplasmic/endoplasmic reticulum calcium ATPase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SERCA1, Calcium pump 1, Calcium-transporting ATPase sarcoplasmic reticulum type, fast twitch skeletal muscle isoform, Endoplasmic reticulum class 1/2 Ca(2+) ATPase
Accession No.	Uniprot ID:O14983
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.