

AT-rich interactive domain-containing protein 4B

Catalog No: #AP80832

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AT-rich interactive domain-containing protein 4B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ARID domain-containing protein 4B, 180 kDa Sin3-associated polypeptide, Sin3-associated polypeptide p180, Breast cancer-associated antigen BRCA1, Histone deacetylase complex subunit SAP180, Retinoblastoma-binding protein 1-like 1, BRCA1, RBBP1L1, RBP1L1, SAP180
Accession No.	Uniprot ID:Q4LE39
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.