

# Aquaporin-2

Catalog No: #AP80838



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Aquaporin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	AQP-2,ADH water channel,Aquaporin-CD,AQP-CD,Collecting duct water channel protein,WCH-CD,Water channel protein for renal collecting duct
Accession No.	Uniprot ID:P34080
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.