

# A disintegrin and metalloproteinase with thrombospondin motifs 2

Catalog No: #AP80841



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADAM-TS 2,Procollagen I N-proteinase,PC I-NP,Procollagen I/II amino propeptide-processing enzyme,Procollagen N-endopeptidase,pNPI,NPI
Accession No.	Uniprot ID:P79331
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.