

Asporin

Catalog No: #AP80878



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Asporin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Periodontal ligament-associated protein 1,PLAP-1
Accession No.	Uniprot ID:Q99MQ4
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.