

Basal cell adhesion molecule

Catalog No: #AP80883



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Basal cell adhesion molecule
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	B-CAM cell surface glycoprotein,Lutheran antigen,CD239
Accession No.	Uniprot ID:Q9MZ08
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.