

Myc box-dependent-interacting protein 1

Catalog No: #AP80895



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myc box-dependent-interacting protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Amphiphysin II,Amphiphysin-like protein,Bridging integrator 1,Amph2,Amphl
Accession No.	Uniprot ID:O08839
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.