

A disintegrin and metalloproteinase with thrombospondin motifs 8

Catalog No: #AP80925



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ADAM-TS 8,METH-2,METH-8,METH2
Accession No.	Uniprot ID:Q9UP79
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.