

A disintegrin and metalloproteinase with thrombospondin motifs 4

Catalog No: #AP80934



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	A disintegrin and metalloproteinase with thrombospondin motifs 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	ADAM-TS 4,ADMP-1,Aggrecanase-1
Accession No.	Uniprot ID:Q9TT93
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.