

A disintegrin and metalloproteinase with thrombospondin motifs 7

Catalog No: #AP80940



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | A disintegrin and metalloproteinase with thrombospondin motifs 7 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | ADAM-TS 7,COMPase |
| Accession No. | Uniprot ID:Q68SA9 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.