

# Cyclic AMP-dependent transcription factor ATF-2

Catalog No: #AP80944



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Cyclic AMP-dependent transcription factor ATF-2                             |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Chicken   |
| Other Names        | cAMP-dependent transcription factor ATF-2,Activating transcription factor 2 |
| Accession No.      | Uniprot ID:O93602   |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                      |

Note: This product is for in vitro research use only and is not intended for use in humans or animals.