

Isoaspartyl peptidase/L-asparaginase

Catalog No: #AP80945



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Isoaspartyl peptidase/L-asparaginase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Asparaginase-like protein 1,Asparaginase-like sperm autoantigen,Beta-aspartyl-peptidase,Glial asparaginase,Isoaspartyl dipeptidase,L-asparagine amidohydrolase,Alp,Gliap,Hiob
Accession No.	Uniprot ID:Q8VI04
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.