

# Myc box-dependent-interacting protein 1

Catalog No: #AP80946



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Myc box-dependent-interacting protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Amphiphysin II,Amphiphysin-like protein,Box-dependent myc-interacting protein 1,Bridging integrator 1,AMPHL
Accession No.	Uniprot ID:O00499
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.